
FRANSOR INDUSTRIE

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CE - 0124

EDITION 2015

DEPARTEMENT DENTAIRE

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Introduction to FransorIndustrie

***FRANSOR INDUSTRIE** is a French metallurgy company set up in 1962 under the name of **FRANCO-SUISSE DENTAIRE S.A.R.L.** It trades precious alloys and materials for dentistry.*

*Since 1969, **FRANCO-SUISSE** has been developing its own range of precious alloys in its 1 800 m² plant near to Paris, which until that time had only been imported from Switzerland.*

Its reputation has grown enabling the company to set up an Industrial Division producing pieces of machinery from precious metals for the Aeronautic, IT, Robotic and Goldsmith industry, using precious metals.

*To reflect its new market positioning in this field, with customers such as **THALES, MATRA, ALCATEL, SAGEM, DASSAULT**, as well its dentistry, activity **FRANCO-SUISSE** was renamed **FRANSOR INDUSTRIE S.A.***

*In March 1994, **FRANSOR INDUSTRIE** was taken over by Brigitte MAIFFRET. She focused on improving the product quality through developing the staff's professional qualifications and the company's customer service.*

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To order gold and precious alloys from **Fransor Industrie**

→ Contact us:

Tél. +33.1.42.42.85.27
Fax +33.1.47.82.89.70
E-mail fransor.industrie@wanadoo.fr
Website: www.fransor-industrie.fr

✱ *We will provide you with **all the** technical information you need*

✱ *We will give you a detailed quote. You will know your exact order value which will be index-linked to the metals current official price on the day of your order.*

→ We deliver within 48 hours !

✱ *We will deliver to your office within 48 hours of receiving your order.*

✱ *Goods will be packaged in a tightly sealed wooden box and sent by post with a declared value for customs.*

✱ *The invoice will be sent with the goods. Payment is due within 2 weeks, by bank transfer.*

✱ *A 2 % discount will be given for immediate payment.*

→ Our Return Policy :

✱ *In the event of a technical problem FRANSOR INDUSTRIE will immediately replace the product. Please contact us immediately to discuss.*

✱ *All our products come with a fourteen day return policy. Please contact us to discuss the issue; as soon as we receive the returned goods, packaged in their original box and undamaged, we can either replace them or issue you with a credit note.*

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FRANSOR INDUSTRIE HAS BEEN AT THE FOREFRONT OF THE DENTAL RESEARCH SINCE 1962

If FRANSOR INDUSTRIE did not exist the world of dentistry would not be what it is today. In the past forty years, FRANSOR INDUSTRIE has been at the forefront of innovation in dentistry, breaking new grounds and developing new innovative products, which have helped dentists and technicians gain valuable time and earn the reputation of being the leading supplier in the French dental industry.

Over the years, FRANSOR INDUSTRIE, has applied its high level of expertise in the dental industry to other industrial activities providing them with precious alloys and precision parts. Thanks to the quality and reliability of its products, today FRANSOR INDUSTRIE is a first choice supplier for all European high tech industries such as the Aeronautic, Aerospatiale, IT, Robotic.

As FRANSOR INDUSTRIE developed its innovation in the industry field, its Dentistry activity has also benefited from these advances in other activities. Over the last few years, FRANSOR INDUSTRIE has focused on developing new solutions for surgery and laboratory, specifically in the technique of attachments, implants and ceramics. This catalogue details the wide range of our dental products. Please do not hesitate to contact us if you require any further information.

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FRANSOR INDUSTRIE **REGISTERED TRADEMARKS**

AZ is the registred trade mark of FRANSOR INDUSTRIE corresponding to the product known under the name of **ASC 52**.

SCOP-STAR is an attachment designed by FRANSOR INDUSTRIE and registered in the Patent Office (INPI) corresponding to the product known under the name of **SCOP**.

FILOPLAST is the registered trade mark of FRANSOR INDUSTRIE corresponding to the product known under the name of **TICOPLAST**.

SG is an attachment designed by FRANSOR INDUSTRIE and registered in the Patent Office (INPI).

MAGNETIZE is an attachment registered by FRANSOR INDUSTRIE in the Patent Office (INPI).

CWJ is an alloy designed by FRANSOR INDUSTRIE and registered in the Patent Office (INPI).

INLOR is an alloy registered by FRANSOR INDUSTRIE in the Patent Office (INPI).

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
ROACH ATTACHMENT...

...A classic to rediscover

- *Two models: 100 % burnable or in steel/burnable*
- *Slotted with correction of looseness for hybrid prosthesis*
- *Absorption of horizontal translation movement and transverse frontal rotation*
- *Roach attachments are sold in boxes of 2.*

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

«AZ»

The “AZ” ball attachment is a small universal joint connector. The indication for the “AZ” (FRANSOR ASC 52) is mainly for precision partial dentures.

- ❖ *The male part is a reversed spring loaded connector allowing adjustable retention.*
- ❖ *Strong and reliable*
NB: no solder
- ❖ *Easy to insert and remove*
- ❖ *Two sizes large and small*


*The AZ exists in two faces: The **super micro** and the **ultra super micro**.*

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
A 100 % burnable female part

The 100 % burnable plastic female part allows inclusion in any fixed prosthesis for better rigidity. It avoids problematic and unreliable solders.

STEP BY STEP

1. After completion of the wax pattern, survey, and trace.
a-the axis of insertion the ridges prepare the abutment teeth to receive the ASC 52
2. With the surveyor and the guiding tool, place the female part on the wax. At this point make sure it is in the proper position regarding both axis and ridges. Place it as close as possible to the abutment teeth (1) and close to the ridge leaving just enough of a gap to avoid gum irritation.
3. Check the lift of the male part (2) there should be a one millimetre gap between both parts on the occlusal.
4. Invest and cast the usual way
5. Before setting up the teeth, fix the male part to the partial frame with one drop of self cure acryl.(3)

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
CONTACTOR STAR

A leader in research, FRANSOR INDUSTRIE created and tested in its laboratories a special plastic which is 100% burnable. It will not leave any residues in the burn out ring. This new material is utilised in all our burn out hinges and attachments, as well as in the new STAR.

❖ *Plastic male part 100% burnable, no residues makes cleaner castings*



❖ *No soldering joint-no messing, no chance to melt the male part*

❖ *Easy to survey, seal with wax on the wax pattern*

❖ *Saves time and money, affordable + easy to do = more savings*

❖ *Cast with your favourite metal, precious or non precious*

❖ *The star is sold complete, the male part can be sold as a separate unit*

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


The first attachment with

a 100% burnable male part

The STAR is a new generation of precision attachment designed for over-dentures. Because of its unique design and small size it can be utilised everywhere, even on the smallest teeth.

DESCRIPTION

- ❖ The male part: is made of a round base with a chamfer ① at the bottom for a better position in the wax pattern, and of a spherical joint in the centre ②. This precision part is made of 100 % burnable plastic: it is also available in precious metal IRIDOR and in non precious metal
- ❖ Female part: is made of a round box with a retentive groove at the top ③ its external superior part includes circular ④ retaining groove, providing a strong fit with the acrylic material. Then 8 active ⑤ leaf spring to ensure a perfect fit of the male part. A plastic ring to avoid stress and wear on the springs and more movement. You can use the included tin washer. The washer will be adapted to the gun and sit between the male part and the female part. The female part is also available in precious metal: IRIDOR, or non precious metal.


STEP-BY-STEP

- 1/Prepare the wax pattern with the post
- 2/Survey the wax pattern to accommodate the STAR in the desire axis
- 3/Place the STAR in the wax pattern and check its proper position and axis. If you use more than one STAR you need to use the STAR survey
- 4/Seal with wax the STAR on the wax pattern
- 5/Invest and cast as usual
- 6/The flasking tool should be used if the pouring or injected acrylic is not done on the master cast
- 7/Final adjustment of the spring is be made with the tool

NB: Do not forget to adjust the springs afterwards

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

CONTACTOR SG

Contactor SG is specially designed for use with removable or hybrid prostheses. It can be sealed in resin with the protective cap, or welded onto a removable prosthesis with the precious alloy threaded cap.

Contactor SG is made of:

- ❖ *A plate with central ball on top.*
- ❖ *A steel or precious metal case which is easily replaceable.*
- ❖ *A sleeve which can be activated through a slot.*
- ❖ *Contactor SG is available in steel, all steel or in precious alloy.*

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


SG and SG Flash

Assembly:

- 1- Position the element in the mouth, and take a mould.
- 2- Pour the mould and mount piece on wax without paying particular attention to the position of the contactor. Then make trial.
- 3- Take plaster key of the vestibular section where the contactor is to be located.
- 4- Suppress the wax and adhere piece onto key.
- 5- Put key on model and mould piece so as to obtain the cavity necessary where the contactor is to be housed with the case.
- 6- Put ball plate in parallel with the rods provided.
- 7- Weld this plate onto flange (never use rod for this operation).
- 8- Insert the tin washer.
- 9- Put brass case on ball with protective cap.
- 10- Fill the empty part with wax and test.
- 11- Loosen brass case, remove sleeve, re-screw brass box.

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
For your over-dentures cases, why not use EZ bar from FRANSOR INDUSTRIE

clip system with an arch-bar you can cast yourself


toproceet of assembly

With FRANSOR'S plastic bar of Ackerman type, it is easy and fast to realise an over-denture. Only one casting and the metal of your choice.
No soldering, it is strong and very simple.

Step-by-step


1. Finish the wax up with the posts in place
2. With a benzene burner gently warm up the plastic bar
3. The bar should be soft and follow the ridge; it should touch every abutment
4. Keep the plastic bar in place till it takes its final shape
5. Attach the bar to each abutment with wax
6. Set up the teeth considering the bar and the ridge
7. Take a matrix of the set up
8. Check the positions of the teeth
9. Place the tin bar over the bar
10. Place the clips over the tin bar
11. The clips should have a light fit
12. Take the teeth set-up to finish;
13. After the flask is open and the wax boiled out slightly cover up the edges of the clips with the stone
14. Finish the acrylic as usual and carefully clean the clips for free movement, adjusting them as necessary.

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
FAST

Burnable-extracoronary-vertical-contactor

The burnable precision attachment FAST is manufactured by FRANSOR INDUSTRIE for precision partial prosthesis. It is especially recommended for three-end saddle types, with one or two saddles.

Step-by-step technique:

- 1) Survey the master cast for undercut ridges, trace the midline and the ridges. Wax-up and survey the sides to receive the male parts of the FAST.
- 2) Using the survey mandrel, place the past male part on the wax pattern, on the axis of the ridge and seal it with wax on it; it should be outside of the crown.
- 3) Invest and cast as usual, with the alloy of your choice.
- 4) Sandblast and clean the casting and finish it at usual.
- 5) The female part should be placed on the male part and the retentive tip of the female part should be fully enclosed in the wax of the saddle set up.
- 6) Re-activate lightly the female part blade if necessary

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Fixed or hybrid prosthesis... **...The solution : 100 % burnout plastic part for attachments**

You have a choice between five 100% burnout models

A Simple type with male part in T-shape. 4.5 mm high.

B The smallest burnout plastic attachment inside the crown. With its stabilisation winglets this model offers a precision connection (non adjustable) 7 mm high.


C Conical intracoronal anchorage, non adjustable. 3.5 mm H

D Simple slotted T-attachment in V-shape with a large retention winglet to include into the wax (6,5 mm high)





E Intracoronal anchorage, non adjustable, with a triangular point. Very good retention in the resin thanks to its perforated winglet. Bridge work with non parallel abutment teeth.

Step by step:

1. Study and survey the master card and establish the path of insertion specific to the patient's mouth.
2. Using wax, attach the male or female part to the wax pattern using the surveyor.
3. Coating will be made with a sable brush then make a slow rise burnout.

For best results and whenever possible, place the female part inside the crown on the distal side on the most mesial abutment tooth.

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

A 100 % burnable all plastic D.E. HINGE



STEP-BY-STEP:


Using wax, stick part A to the pre-wax pattern without over contouring on either side, respecting the axis of the ridge.

1. Position the D.E. hinge on the model
 2. Stick the A part on the wax pattern
 3. Remove the axle from the B part
 4. Cast the B part separately
 5. After casting, adjust the B part on the A part
- Insert the axle in the D.E. hinge and process with the teeth set-up as usual

Note: After casting sandblast and finish both parts on the master cast. After the teeth are set and the acrylic cured, carefully warm up the metal parts to free the movements of the D.E. hinge.

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



SCREW for removable or irremovable core frame

Taper thread screw in Iridor 1480°C with collar and cylindrical tube. It's recommended for any application on semi fixed bridges in platinum-coated gold or auto ceramic.

Assembly:

1. *Insert the lubricated sleeve into the wax pattern of the element(s). Stop when it's at the same level with the sleeve.*
2. *Remove the sleeve and replace it by a steel shank.*
3. *Apply the coating and cast*
4. *Once the casting is ready, remove the steel shank out.*
5. *Bevel the orifice off $\frac{1}{2}$ mm with a round chamfer cut the diameter of which must be larger than $\frac{1}{2}$ mm.*
6. *Insert the sleeve by sliding frictions*



IMPORTANT

7. *Before soldering fill the sleeve with graphited paste.*
8. *Put 2 or 3 small pieces of soldering into the milling and solder*
9. *Clean inside of the sleeve with ultra-son.*
10. *Insert the screw into the hoop, screw it home then half a turn before shaping the removable part. Lubricate the casted part embodying the sleeve and shape with the amorphous wax.*
11. *Take the screw out and cast the shaped part.*
12. *Finish the removable element.*

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grandeur nature

MAGNETIZE

**THE MAGNETIC ATTACHMENT
OF FRANSOR INDSTRIE (CE 0124)
IS ON SALE AT YOUR RETAILER**

The **MAGNETIZE** kit includes:


- ☞ 1 magnet FRANSOR INDUSTRIE
- ☞ 1 washer HOSTOPLAST
- ☞ 1 resilient washer

Thanks to its, **MAGNETIZE** is used in
Most of clinical cases requiring a magnet
attachment.

TECHNICAL CHARACTERISTICS

- ☞ Height: 2.25 mm
- ☞ Diameter: 4.60 mm
- ☞ Shape: 300 gr
- ☞ Material : Inox 316 L
(for medical quality)

ITS ASSETS !

The **HOSTOPLAST** washer once fitted
on the wax enables the technician to
make a very flat plateau. Its mechanical
properties allow a total anchorage in
the wax.

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The famous cylindrical-conical posts

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Calcination, without any residue.

- The additional advantage to of rigidity and flexibility.

-Useful in surgery as well as in laboratory:

IN SURGERY:

when the standardized posts are not suitable for an impression

IN LABORATORY:

To cast personalized posts or shanks. Their 8 cm length resolves any problem.

PACKAGING:

Each tube takes shanks which have the same diameter.

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
FRANSOR INDUSTRIE

The first ones to deserve the label « high tech » !

Conceived in close collaboration with dental professionals, our dental gold are melted in our manufactory division near to Paris, a genuine reference of the high-tech industry.

The main properties of our gold production are the alloy purity and a range of products really suited to dental care and today's techniques. We assure competitive prices and prompt service.

Our laboratory maintains tight control on the products which do conform to the most stringent norms.

FRANSOR INDUSTRIE is the first independent French supplier in the technical gold market. Thanks to our highly qualified staff and its investment in the latest technologies, we are able to provide the best quality-price ratio.

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
Cwj

Platinum and Rhodium are introduced into the composition to increase resistance to abrasion, dulling and corrosion by ten fold. Without silver nor Palladium and mechanical characteristics suits all the ceramics available on the market.

Ceraminor

Gold alloy, white in colour, without Silver. It has the highest young's modules and mechanical qualities, compatible with the porcelain available on the market. It suits drilling and affiliated works. Very easy with primary soldering.

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inlaycast

Because of its hardness, this lovely yellow alloy is ideal for inlays.

It is very rich in precious metals and safe
It is self-hardening and is only produced in ingots.

crownor C

It is a pinkish alloy which contains significant amount of platinum and palladium. It is self-hardening, very homogeneous and fine-grained. It is supplied in the form of ingots and wires of all sizes.

It is recommended for framework prostheses which have to withstand great masticatory stress, as well as crowns, bridges and clasps.

platinor C

A pinkish yellow self-hardening alloy containing platinum and metals of the platinum group.

It is 15 % lighter than the classical alloys and is very strong.

It is produced in the form of ingots, and is specially designed for bridges, crowns and framework prostheses.

Seulor C

Pale yellow alloy, auto hardening, elastic. It is a Palladium gold alloy suitable for the requirements of today's market.

It is one of the most reliable alloys, 17% lighter than classical alloys and 6 times stronger. It is ideally suite for thin pieces.

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Palgor C

White grey Alloy, self-hardening, modern, with a higt precious metals content. The most sought after and appreciated economical white metal on the market.

Palgor provides resilience to change and resistance. In addition due to its lightness, allows the manufacture of economical prosthesis such as inlays, bridges, crowns etc.

Castor C

White grey Alloy, very economical, 30% lighter than 22 carat gold. Self-hardening, it allows the manufacture of mocile or fixed prothesis. It replaces steel advantageously and is easy to manipulate. It is recommended to cast the cylinder quickly and not go over 200C.

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